

Splash screen

Content

1. Dia-Buddy (text)
-

On boarding 1 screen

Content Elements

1. Woman Eating (Image)
2. Track your sugar levels (Text)
3. Progress indicator (shape)

Function Elements

1. Next button (button)
-

On boarding 2 screen

Content Elements

1. People Supporting (Image)
2. Join a diabetes community (Text)
3. Progress indicator (shape)

Function Elements

1. Next button (button)
-

On boarding 3 screen

Content Elements

1. Balanced meal (Image)
2. Monitor meals & Carbs (Text)
3. Progress indicator (shape)

Function Elements

1. Get Started (button)

Sign up screen

Content Elements

1. Sign Up (text)
2. Create your account to start tracking your health (text)
3. Or Sign up with (text)

Function Elements

1. Back button (icon)
 2. Name (field)
 3. Email (field)
 4. Password (field)
 5. Confirm password (field)
 6. Agree to terms & conditions (checkbox)
 7. Continue (button)
 8. Have an account? Login (link)
 9. Google (icon)
 10. Facebook (icon)
 11. Apple (icon)
-

Phone number screen

Content Elements

1. Phone Number (text)
2. Enter your phone number to verify your account (text)
3. Progress indicator (circular shape)

Function Elements

1. Back button (button)
 2. Phone number (field)
 3. Next button (button)
-

Account type screen

Content Elements

1. Select account type (text)
2. Select whether you are a patient or a caregiver (text)
3. Patient card (Image + label)
4. Care giver card (Image + label)
5. Progress indicator (circular shape)

Function Elements

1. Back button (button)
 2. Card selection (interactive cards)
 3. Continue button (button)
-

Caregiver type screen

Content Elements

1. Select caregiver role (text)
2. Select whether you are a doctor or a family member (text)
3. I am a Doctor Card (image + label)
4. I am a Nurse/Family member Card (image + label)
5. Progress indicator (circular shape)

Function Elements

1. Back button (button)
 2. Card Selection (interactive cards)
 3. Continue button (button)
-

Gender screen

Content Elements

1. Gender (text)
2. This will help us personalize your experience (text)
3. Female image (image)
4. Male image (image)
5. Progress indicator (circular shape)

Function Elements

1. Back button (button)
 2. Gender selection (image)
 3. Continue button (button)
-

Height screen

Content Elements

1. Height (text)
2. Choose your preferred unit and provide your height (text)
3. Progress indicator (circular shape)
4. Ruler icon (icon)

Function Elements

1. Back button (button)
 2. Cm (button)
 3. Ft (button)
 4. Height input (ruler slider)
 5. Continue button (button)
-

Weight screen

Content Elements

1. Weight(text)
2. Choose your preferred unit and provide your weight (text)
3. Progress indicator (circular shape)
4. Weighing scale (icon)

Function Elements

1. Back button (button)
 2. Kgs (button)
 3. Lbs (button)
 4. Weight input (circular slider)
 5. Continue button (button)
-

Birth date screen

Content Elements

1. Birthdate (text)
2. This helps us personalize your experience (text)
3. Birthdate (icon)
4. Progress indicator (circular shape)

Function Elements

1. Back button (button)
 2. Birthdate selector (date picker)
 3. Continue button (button)
-

Blood type screen

Content Elements

1. Blood Group (text)
2. Select your blood group (text)
3. Blood drop (icon)
4. Progress indicator (circular shape)

Function Elements

1. Back button (button)
 2. A+ (button)
 3. A- (button)
 4. B+ (button)
 5. B- (button)
 6. O+ (button)
 7. O- (button)
 8. AB+ (button)
 9. AB- (button)
 10. Continue button (button)
-

Diabetes type screen

Content Elements

1. Diabetes type (text)
2. What kind of diabetes do you have (text)
3. Progress indicator (circular shape)

Function Elements

1. Back button (button)
 2. I have Type 1 diabetes (button)
 3. I have Type 2 diabetes. (button)
 4. I have gestational diabetes. (button)
 5. I'm not sure about my type. (button)
 6. Other (field)
 7. Continue button (button)
-

Medicine type screen

Content Elements

1. Medications (text)
2. What type of medicine do you take (text)
3. Syringe card (label + icon)
4. Pills card (label + icon)

Function Elements

1. Back button (button)

2. Medication type selection (interactive cards)
 3. Medicine search (field)
 4. Frequency intake (drop down)
 5. Sign Up (button)
-

Log in screen

Content Elements

1. Log in (text)
2. Access your account to continue your journey (text)
3. Or log in with (text)

Function Elements

1. Email (field)
 2. password(field)
 3. Keep me logged in (checkbox)
 4. Forget password (link)
 5. Log in (button)
 6. Google (icon)
 7. Facebook (icon)
 8. Apple (icon)
 9. DO not have an account? Sign up (link)
-

Scan Qr Code screen

Content Elements

1. Connect to your Caregiver(text)

Function Elements

1. Switch between scan and generate Qr code (Tab selector)
 2. Qr Scanning Area (camera input)
 3. Connect button (button)
 4. Skip (link)
-

Generate Qr Code screen

Content Elements

1. Connect to your Caregiver(text)
2. Qr Code display (image)

Function Elements

1. Switch between scan and generate Qr code (Tab selector)
 2. Connect button (button)
 3. Skip (link)
-

Streak screen

Content Elements

1. Streak fire (image)
2. Ongoing streak Day (text)
3. How are you feeling today? (text)

Function Elements

1. Sad face (icon)
 2. Moderate face (icon)
 3. Happy face (icon)
-

Dashboard screen

Content Elements

1. Greeting (text)
2. Glucose level Diagram card (Label + Diagram + image)
3. Carbs level Diagram card (Label + Diagram + image)
4. Heart Rate Card (Card + Icon + Value + Label)
5. Upcoming Medications
6. Morning (icon)
7. Morning (text)
8. Evening(icon)
9. Evening (text)

Function Elements

1. Notification bell (icon)
2. Streak (icon)
3. Glucose or carbs diagram selection (interactive cards)
4. Carbs and glucose and heart rate (Interactive Cards)
5. Ai chatbot (icon)
6. Emergency (button)
7. See more (link)
8. Medicine (checkbox)
9. Navigation Bar
 - A .home (icon)
 - B. analytics (icon)
 - C. Plus (icon)
 - D. community (icon)
 - E. profile (icon)

Blood glucose entry screen

Content Elements

1. Glucose (text)

Function Elements

1. Save (button)
2. Cancel (button)
3. Glucose level (field)
4. Blood pressure (field)
5. Date (picker)
6. Time (picker)
7. Measuring time (picker)

Food entry screen

Content Elements

1. Food (text)
2. Choose saved meal (text)

3. Card meal (image+ label)

Function Elements

1. Save (button)
 2. Cancel (button)
 3. Meal Type (dropdown)
 4. Meal Name (field)
 5. Date (picker)
 6. Time (picker)
 7. Photo scan (link)
 8. Units (radio buttons)
 9. Saved meal cards (Interactive cards)
-

Medications screen

Content Elements

1. Medications (text)
2. Notification bell (icon)
3. Medication progress today (circular shape)
4. Morning (icon)
5. Morning (text)
6. Evening(icon)
7. evening(text)
8. Medication card (icon + label + time)

Function Elements

1. Back button (button)
 2. Notification (toggle)
 3. Medication card (checkbox)
 4. Plus icon (icon)
-

Add medication screen

Content Elements

1. Rounded rectangle (shape)
2. Medication (text)
3. Time (text)

4. Medication type (icon)

Function Elements

1. Edit (button)
 2. Delete (button)
 3. Medication box (checkbox)
-

Emergency Click screen

Content Elements

1. Location of caregiver (Maps)
2. The Caregiver is on their way (text)
3. Caregiver name (text)
4. Time remaining(text)
5. Distance in Km (text)

Function Elements

1. Back button(button)
 2. Ping (button)
-

Analytics day screen

Content Elements

1. Analytics (text)
2. Day and date (text)
3. Mood (icon)
4. Carb and glucose entries(graph)
5. glucose levels (graph)
6. crab levels (graph)

Function Elements

1. Back button(button)
2. Food entries (icon)
3. Food entry details (interactive card)
4. Glucose level details (interactive card)
5. day, week, month(tabs)
6. edit mood (icon)
7. Navigation Bar
 - A .home (icon)
 - B. analytics (icon)
 - C. Plus (icon)
 - D. community (icon)
 - E. profile (icon)

Analytics week screen

Content Elements

1. Analytics (text)
2. Mood tracker (graph)
3. Glucose level tracker (graph)
4. Carb level tracker (graph)
5. week date range (text)

Function Elements

1. back button(button)
2. Shuffle between weeks (icon)
3. Navigation Bar
 - A .home (icon)
 - B. analytics (icon)
 - C. Plus (icon)
 - D. community (icon)
 - E. profile (icon)

Analytics month screen

Content Elements

1. Month name (text)
2. Analytics (text)

Function Elements

1. Back button(button)
 2. day, week, month (tabs)
 3. Days display (calendar)
 4. Navigation Bar
 - A .home (icon)
 - B. analytics (icon)
 - C. Plus (icon)
 - D. community (icon)
 - E. profile (icon)
-

Community screen

Content Elements

1. Posts (text)
2. Recent chats (text)

Function Elements

1. see all (link)
2. Post cards (interactive cards)
3. Add post (icon)
4. Chats (icon)
5. see all (link)
6. I Navigation Bar
 - A .home (icon)
 - B. analytics (icon)
 - C. Plus (icon)
 - D. community (icon)
 - E. profile (icon)

All Posts

content elements

1. All posts (text)

function elements

1. post cards (interactive cards)
2. back button (button)
3. new post (button)
4. drafts (button)

Create posts

content elements

1. New post (text)
2. Account name (text)
3. account icon (image)

function elements

1. cancel (icon)
2. attachment (icon)
3. send (icon)
4. Save as draft (button)
5. Post text (field)

Drafts

content elements

1. Drafts(text)

function elements

1. back button (button)
 2. post cards (interactive cards)
-

Community

content elements

1. chat (text)
2. My communities (text)

Function elements

1. back button (button)
 2. Community , Drs (tabs)
 3. Search (field)
 4. New community (icon)
 5. chats (icon)
-

Drs

content elements

3. chat (text)
4. My Drs (text)

Function elements

6. back button (button)
 7. Community , Drs (tabs)
 8. Search (field)
 9. Add dr (icon)
 - 10.chats (icon)
-

On Dr Profile click

Content elements

1. Dr profile (text)

2. Profile picture (image)
3. Dr name (text)
4. Availability (text)
5. experience years (text)
6. rating (text)
7. Qualifications (text)
8. clinics (text)

function elements

1. back button (button)
2. Message (button)

Streak Calendar

content element

1. Streak Calendar (text)
2. Streak (text)
3. Fire (image)
4. Details (text)
5. Month (text)

function element

1. back button (button)
2. Notifications (icon)
3. month shuffle (icon)
4. days display (calendar)

Patient profile screen

Content Elements

1. Profile (text)
2. Profile card (name + email + profile image)
3. Glucose type (card + label)
4. Blood type (card + label)
5. Height (card + value)

6. Gender (card + value)
7. Weight (card + value)
8. Age (card + value)

Function Elements

1. Edit icon (icon)
 2. Back button (button)
 3. Log out (icon)
 4. Settings (card)
 5. Caregiver (card)
 6. Medicine (card)
 7. Device connected (card)
-

Update Profile screen

Content Elements

1. Update profile (text)

Function Elements

1. Back button (button)
 2. Change profile image (image)
 3. Change email (field)
 4. Change name (field)
 5. Change glucose type (field)
 6. Change height (field)
 7. Change weight (field)
 8. Change birthdate (field)
 9. Change blood type (field)
 10. Change medicine (field)
 11. Add new caregiver (plus icon)
 12. Add another caregiver (plus icon)
 13. Remove caregiver (minus icon)
 14. Save changes (button)
 15. Discard (button)
-

Successful update profile screen

Content Elements

1. Shape (animation)
 2. Data changes successfully (text)
-

Devices screen

Content Elements

1. My devices (text)
2. My connected devices (text)
3. Other devices (text)

Function Elements

1. Available to connect (cards)
 2. My devices (cards)
 3. Add new device (button)
 4. Back button (button)
-

Connecting device screen

Content Elements

1. Connecting your device (text)
2. Please keep your device close and wait while we set things up (text)
3. Loading shape (animation)
4. This may take a few seconds (text)
5. My device (image)

Function Elements

1. Connect button (button)
 2. Back button (button)
-

Notifications screen

Content Elements

1. Notifications (text)
2. Bill (icon)

Function Elements

1. All (tab)
2. Community (tab)
3. Medications (tab)
4. Doctors (tab)
5. Notification (card)
6. Delete (button)
7. Back (button)

Help center screen

Content Elements

1. Help center (text)

Function Elements

1. Back button (button)
2. Chatbot (list item)
3. FAQ (list item)
4. Info & Advice (list item)
5. Terms of services (list item)

FAQ screen

Content Elements

1. FAQ (text)
2. Find quick answers to common questions? (text)

Function Elements

1. Card question
 - A. How do I connect my glucose monitoring device? (Q text)
 - B. Arrow (icon)
 - C. Go to Devices > Add New Device and select your device from the list. Follow the on-screen instructions to pair it via Bluetooth (A text)
2. Card question
 - A. My device is not syncing. What should I do? (Q text)
 - B. Arrow (icon)
 - C. Make sure Bluetooth is turned on and your device is charged. If the problem continues, restart both your phone and the device, then try again (A text)
3. Card question
 - A. How do I log my glucose manually? (Q text)
 - B. Arrow (icon)
 - C. Tap the '+' icon on the navigation bar, select Glucose Entry, and type your reading (A text)
4. Card question
 - A. Can I track my daily calorie intake? (Q text)
 - B. Arrow (icon)
 - C. Yes. Use the '+' button to add meals and the app will automatically calculate calories and carbs for you (A text)
5. Card question
 - A. Is there a support community I can join? (Q text)
 - B. Arrow (icon)
 - C. Yes, you can access the Community tab to connect with other users, share experiences, and get support (A text)

6. Card question

A. How do I report inappropriate behavior in the community? (Q text)

B. Arrow (icon)

C. Open the user's profile, tap the 'Report' option, and select the reason. Our team will review it quickly (A text)

7. Card question

A. How can I edit my profile information? (Q text)

B. Arrow (icon)

C. Go to Profile > edit icon and select the field you want to update (name, age, gender, etc.) (A text)

8. Card question

A. I forgot my password. How can I reset it? (Q text)

B. Arrow (icon)

C. On the login screen, tap Forgot Password and follow the steps to reset it via email (A text)

9. Card question

A. How do I add my caregiver or doctor to view my data? (Q text)

B. Arrow (icon)

C. Go to Settings > Caregiver Access and generate a QR code. Your caregiver can scan the code using their app to connect instantly (A text)

10. Card question

A. Can I remove a caregiver after adding them? (Q text)

B. Arrow (icon)

C. Yes. In Settings > Caregiver Access, select the caregiver and tap Remove Access. This will immediately revoke their connection (A text)

11. Back button (button)
12. Search field (field)

Terms of services screen

Content Elements

1. Terms of services (text)
2. User Responsibilities (text)
 - a. You are responsible for entering accurate health data and keeping your login credentials secure (text)
3. Data Privacy & Security (text)
 - a. Your personal and health information is encrypted and never shared without your consent (text)
4. Data Privacy & Security (text)
 - b. Your personal and health information is encrypted and never shared without your consent
5. Community Guidelines (text)
 - a. Respect others in the support community. No spam, offensive, or harmful content is allowed (text)
6. Limitations & Disclaimers (text)
 - a. The app provides tracking and educational tools but does not replace professional medical advice (text)
7. Contact Information (text)
 - a. For support, contact us at support@diabuddy.com or visit Help Center within the app (text)

Function Elements

1. Back button (button)
 2. support@diabuddy.com (link)
-

Info & Advice screen

Content Elements

1. Info & Advice (text)
2. Monitor Blood Sugar Regularly (text)
 - a. Check your glucose levels daily using your device to understand how food, exercise, and medication affect your body (text)
 - b. Blood test (image)
3. Eat Balanced Meal (text)
 - a. Focus on vegetables, lean protein, whole grains, and avoid high-sugar foods to maintain stable glucose levels (text)
 - c. Healthy meal (image)
4. Stay Active (text)
 - a. Regular exercise, even a 20-minute walk, helps your body use insulin more effectively (text)
 - a. Fitness activity (image)
5. Stay Hydrated (text)
 - a. Drink enough water throughout the day to prevent dehydration and help regulate blood sugar levels (text)
 - b. Water glass (image)
6. Take Medications on Time (text)
 - a. Always follow your doctor's instructions and set reminders to take your medicines or insulin on schedule (text)
 - b. Medicine Reminder (image)
7. Manage Stress (text)
 - a. Stress can affect blood sugar. Try relaxation exercises like deep breathing, yoga, or meditation (text)
 - b. Meditation (image)
8. Get Regular Checkups (text)
 - a. Schedule routine checkups with your doctor to track your progress and adjust your treatment plan if needed (text)
 - b. Doctor visit (image)

9. Join a Support Community (text)

- a. Share experiences with others managing diabetes. Support can motivate and help you stay on track (text)
- b. Community support (image)

Function Elements

1. Back button (button)
-

AI chatbot screen

Content Elements

1. DiaBot (text)
2. Let me help you today? (text)
3. AI bot (image)

Function Elements

1. Back button (button)
 2. Ask (button)
-

AI chat screen

Content Elements

1. DiaBot (text)
2. DiaBot avatar (image)
3. User question (text)
4. Bot answer (text)

Function Elements

1. Back button (button)

2. User input field (text field)
3. Send button (button)
4. Quick suggestions (cards)

Caregiver management screen

Content Elements

1. My caregivers (text)
2. Caregiver card (card)

Function Elements

1. Back button (button)
2. Add new one (button)
3. Delete icon (icon)

Add caregiver screen

Content Elements

1. Add caregivers (text)

Function Elements

1. Back button (button)
2. Add caregiver (button)

Forgot password screen

Content Elements

1. Forgot password (text)
2. Enter your registered email to reset your password (text)
3. Already have an account (text)

Function Elements

1. Email (field)
 2. Send code (button)
 3. Log in (link)
-

Verification screen

Content Elements

1. Check your email (text)
2. We've sent a password reset link to your email address. Please check your inbox or spam folder (text)

Function Elements

1. Open email app (button)
 2. Didn't get the email please Resend (link)
-

Reset password screen

Content Elements

1. Reset your password (text)
2. Enter a new password (text)
3. Password Strength Indicator (progress bar)

Function Elements

1. New password (field)
 2. Confirm password (field)
 3. Save new password (button)
-

Success Confirmation screen

Content Elements

1. Shape (animation)
2. Password reset successful (text)

Function Elements

1. Log in (button)